

ReTurn7400™

SystemRoMedic™

Patented design
Low rising ladder,
ideal for shorter adults
and children
Activates the user
Eliminates heavy lifting
for the caregiver*



Safe and active sit-to-stand and transfer – without any heavy lifting

ReTurn7400 is an easy-to-use and incredibly effective assistive device for sit-to-stand and transfer, but it is also a rehabilitation device. ReTurn7400 enables safe and active sit-to-stand and transfer to or from bed, wheelchair or toilet. ReTurn7400 can also be used to facilitate repositioning farther back in a chair or wheelchair. For the caregiver, ReTurn7400 enables easy transfers without any heavy lifting. For users who still retain some degree of strength, ReTurn7400 stimulates the natural pattern of movement during sit-to-stand, while at the same time strengthening the muscles and functional ability. This makes it the ideal assistive device in a rehab unit, and a user who uses ReTurn7400 at home can normally live at home longer with the help of only one caregiver.



handicare

ReTurn7400



ReTurn7400 and ReTurn7500 have the same stable base plate, which makes it possible to change between low and standard-height rising ladders, and thereby cover a broader range of needs and get the best use out of the product.

Product specification

Item No.	Product	Description
7400	ReTurn7400	Transfer platform, complete with low rising ladder
7401	Rising ladder, low	For ReTurn7400 with serial numbers from 373317697400500001 and for ReTurn7500 with serial numbers from 373317697500500001
7501	Rising ladder, standard height	For ReTurn7500 with serial numbers from 373317697500500001 and for ReTurn7400 with serial numbers from 373317697400500001
7501i	Rising ladder, standard height	With an opening for hooking of ReTurnBelt, for ReTurn7500 with serial numbers from 373317697500500001 and or ReTurn7400 with serial numbers from 373317697400500001



Different models of ReTurn provide solutions for all situations and users

The basic principle behind ReTurn is that the user's own strength and ability must be utilized and improved. All ReTurn models therefore help to activate the user and encourage the user's natural patterns of movement.



More information about the different models of ReTurn and all accessories can be found at www.handicare.com/us.

Easy-to-use and very versatile

ReTurn7400 is a manual sit-to-stand and transfer aid that is very easy to understand and use, both for the user and the caregiver. ReTurn7400 is also a versatile product that can be used in several different ways to safely and efficiently assist a user during sit-to-stand and positioning in a seated position, as well as for short transfers indoors. In addition, ReTurn7400 is also an excellent training device, making it an ideal aid to the rehabilitation process.

ReTurn7400 has a low rising ladder that is ideal for shorter adults and children. It also works well in special care situations, such as in an x-ray unit. ReTurn7400 is intended for users weighing max. 150 kg/330 lbs.

Incredibly easy to manoeuvre

ReTurn7400 has a unique, patented* design which makes the transfer platform incredibly easy to manoeuvre. ReTurn7400 rolls, turns and pivots smoothly and easily on all surfaces, even in confined spaces. The device is easily positioned perfectly, close to a bed, wheelchair or toilet, and moving ReTurn7400 requires minimal exertion, even with a heavier user.

Stable and secure base plate

A low centre of gravity and a broad underlying support surface, with wheels at all corners, make the base plate stable and safe. Two large, fixed and lockable centre wheels make ReTurn7400 easy to roll and directionally stable, while the castors at the corners make it very easy to manoeuvre.

The rising ladder allows space for a natural pattern of movement and closeness

ReTurn7400 is equipped with a rising ladder that is fully open at the bottom, while the upper section allows many possibilities for grasping the ladder, including a central 'power grip'. The rising ladder leans forward, away from the user, allowing the user plenty of space to get really close to the assistive device, so as to be able to 'climb' and/or pull himself up while rising to a standing position. The user is also supported by the ladder during standing training and transfer. The design allows plenty of work space also for the caregiver, who can remain very close and help to support the user from different directions.

The lower leg supports can be adjusted both vertically and laterally. They can also be angled, making it easy to individually adapt the settings according to the user's height and other factors. This function also makes it possible for several individuals to use the same ReTurn7400.

ReTurn7400 is easy to assemble and disassemble for transport and storage.

Technical information

SWL	150 kg / 330 lbs
Overall height	995 mm / 39.17"
Base plate, length x width x height	685 mm x 570 mm x 40 mm / 26.97" x 22.44" x 1.57"
Weight	16.5 kg / 36.38 lbs

* Patent: SE-0300528-7 • Registered Community Design: RCD-529631-1

Accessories

By combining ReTurn7400 with different accessories the user can gain additional support and a greater sense of security during sit-to-stand, standing training and transfer. The choice of accessories depends on the situation and the user's needs and functional ability.

Item No.	Product	Description
6031--6038	ReTurnBelt**	Size XXS--XXL
6133--6136	ReTurnBelt, wipeable	Size S--XL
6233--6236	ReTurnBelt, disposable, non-woven	Size S--XL
6061	SupportStraps	Accessory for belts
7250	HeelStrap	For ReTurn7400/7500
7215	Padding for leg support	For ReTurn7400/7500
0954	Bag	For ReTurn7400/7500

** Patent: SE-0401517-8



handicare

2201 Hangar Place, Suite 200
Allentown, PA 18109
Tel: 610 266 5260 • Fax: 610 266 5266
Toll Free: 866 276 LIFT/5438
www.handicare.com/us
customerservice.usa@handicare.com

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

Transfer, between two locations. **Positioning**, within the same location.

Support, for mobility. **Lifting**, both manual and mechanical.