# EasyGlide

SystemRoMedic™

User manual - English

Bruksanvisning - Svenska

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**REF**

5050

5020

5060

5070

**SWL: 135 - 150 kg/297 - 330 lbs**

**SWL: 135 - 170 kg/297 - 375 lbs (5050)**

## Table of contents

<b>English</b> .....	3-4
<b>Svenska</b> .....	5-6
<b>Norsk</b> .....	7-8
<b>Dansk</b> .....	9-10
<b>Soumi</b> .....	11-12
<b>Deutsch</b> .....	13-14
<b>Nederlands</b> .....	15-16
<b>Français</b> .....	17-18
<b>Italiano</b> .....	19-20
<b>Español</b> .....	21-22
<b>Size</b> .....	23

The sliding board EasyGlide is an ingenious invention. The board helps transfer sitting users to and from wheelchairs, beds, toilets, car seats and many other seated positions. It is especially suitable for levelling out height differences and for bridging short distances.

The rounded corners and the smooth, low-friction top surface provide for easy sliding transfers and make the sliding board very easy to position and remove. The double anti-slip tape on the underside prevents the board from sliding against the underlying surface during transfer.


## Check Safety

### Visual inspection

Check the condition and function of the assistive device regularly. Verify that the materials are not worn or faded. If there are signs of wear, the product must be discarded.

### Stability test

Carefully bend the EasyGlide to check strength and durability. The board should be both elastic and rigid. Discard the EasyGlide when its strength and durability are diminished.

## Always read the user manual

Always read the manuals for all assistive devices used during a transfer.

Keep the manual where it is accessible to users of the product.




Always make sure that you have the right version of the manual.

The most recent editions of manuals are available for downloading from our website, [www.handicare.com](http://www.handicare.com).


## Important information

Do not leave the user unattended during a transfer situation.

## Placement/removal of EasyGlide


Place the wheelchair at an angle towards the side of the bed, and as close to the bed as possible. Remove the armrest that is closest to the bed and remove leg supports. Place the EasyGlide like a bridge between the bed and the wheelchair. If it is an EasyGlide model on which one long side can be bent, this should be placed against the wheelchair wheel.

To position the EasyGlide under the patient, lean the patient slightly in the opposite direction.

When the transfer is complete, EasyGlide is removed by bending upwards.

## Using EasyGlide

Active users can transfer themselves with the help of EasyGlide. Grasp the armrest and lean over towards the wheelchair. It is important to remember to lean once more towards the bed when the transfer is complete, so as to be able to move back in the middle of the chair and slide all the way back.

If the patient requires support and assistance during the transfer, the EasyGlide should be complemented with an EasyBelt/FlexiBelt.

The caregiver can work from behind or in front of the patient.

Tips! Our FootStool can be used to advantage to raise the knees

EasyGlide can also be used for transfers into and out of a car.

## Material

Polyethen

## Care of the product

The surface of the product can be cleaned with warm water or a soap solution. If disinfection is needed, 70 % ethanol, 45 % isopropanol or similar should be used.

**CE** Medical Device Class I

The product complies with the requirements of the Medical Devices Directive 93/42/EEC.

Contact your local distributor if you have any questions about the product and its use. See [www.handicare.com](http://www.handicare.com) for a complete list of distributors.

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## Simple solutions for great results

SystemRoMedic™ is the name of Handicare's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic™ is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

## Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic™ is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic™ offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website; [www.handicare.com](http://www.handicare.com).

## For questions about the product and its use

Please contact your local Handicare and SystemRoMedic™ representative. A complete list of all our partners with their contact details can be found on our website; [www.handicare.com](http://www.handicare.com).

### ***Handicare offers solutions and support to increase the independence of disabled or elderly people as well as to improve the convenience of those who are caring for them.***

*The Handicare Group is one of the leading healthcare companies in Europe with its own manufacturing organizations and sales companies in Norway, Sweden, Denmark, Germany, the Netherlands, Great Britain, France, China, Canada and the USA. Handicare's products are also distributed by partners in more than 40 countries worldwide. Our wide range of high-quality products includes various manual and power wheelchairs, seating systems, scooters, a complete easy transfer system and other patient handling aids, stairlifts, car adaptations, rise and recline chairs, and bathing and toileting products.*

**handicare**

Handicare AB  
Maskinvägen 17  
SE-972 54 Luleå, SWEDEN  
Tel: +46 (0)8-557 62 200  
Fax: +46 (0)8-557 62 299  
E-mail: [info@handicare.se](mailto:info@handicare.se)  
[www.handicare.com](http://www.handicare.com)