


Lifting slings guide

All slings = SWL 300 kg



Model/support	Low back +	Shoulder high ++	High back +++ / ++++		
Opening/ leg supports Large opening divided leg supports	HygieneLBSling Polyester	HygieneSling Polyester, with belt Polyester, without belt Polyester, with belt and extra back support	HygieneHBSling Polyester		
 Standard opening divided leg supports	LowBackSling Polyester Polyester net	BasicSling Polyester Polyester net	HighBackSling Polyester Polyester net Plastic net Disposable, non-woven	ShellHBSling Polyester net	FlexibleSling Polyester net
 Small opening divided leg supports		AmpSling Polyester Polyester net	AmpHBSling Polyester		
 No opening Undivided leg supports with hygiene opening		ClassicSling Polyester net	ClassicHBSling Polyester net	FlexibleSling, undivided Polyester net	FlexibleSling, undivided with hygiene opening Polyester net

● ● ● ● ● ● ● Handicare lifting slings are available in sizes ranging from XXS to XXL. The smallest sizes fit children and small adults, and the largest sizes fit large and plus size users. The three sizes in the middle of the scale, S, M and L, fit most normal-sized adults. We also have special lifting sling models in XXXL and XXXXL for extremely large and heavy users. A more detailed overview, with article numbers and measurements for all sling models in all materials, is featured on our website.

X X X S M L X X X
 X X S L X X
 X S L X
 S L