

Technical drawings

Shower chair textile seat, mobile


Article no:

LI2140.000 - Shower chair textile seat, mobile without footrests, 540 mm

LI2140.100 - Shower chair textile seat, mobile without footrests, 490 mm

LI2141.000 - Shower chair textile seat, mobile, 540 mm

LI2141.100 - Shower chair textile seat, mobile, 490 mm



	A	B
LI2140.000/ LI2141.000	540 mm/21.3"	520 mm/20.5"
LI2140.100/ LI2141.100	490 mm/19.3"	470 mm/18.5"

	C
LI2141.000	280-410 mm/11"-16.1"
LI2141.100	280-360 mm/11"-14.2"

	mm	inch
D	900	35.4"
E	450	17.7"
F	490	19.3"
G	415	16.3"
H	250	9.8"
I	625	24.6"
J	450	17.7"