

Technical drawings


Shower chair textile seat, self-propelled

Article no:

LI2142.000 - Shower chair textile seat, self-propelled, 540 mm

LI2142.100 - Shower chair textile seat, self-propelled, 490 mm

LI2142.108 - Swimming pool hoist chair, self-propelled



	A	B	C
LI2142.000	540 mm/21.3"	520 mm/20.5"	280-410 mm/11"-16.1"
LI2142.100/LI2142.108	490 mm/19.3"	470 mm/18.5"	280-360 mm/11"-14.2"

	mm	inch
D	1080	42.5"
E	450	17.7"
F	490	19.3"
G	415	16.3"
H	250	9.8"
I	650	25.6"
J	450	17.7"