

## For sit-to-stand and transfers with turning.

TurnSafe is an ergonomic and easy-to-use aid for sit-to-stand and transfers with turning. TurnSafe activates and supports users with impaired balance and strength in the legs during transfers from bed to wheelchair or from wheelchair to the toilet. The turning plate is small and convenient and very easy to manoeuvre. It's an excellent aid which spares the backs of the carers and reduces the risk of work related injuries.



TurnSafe is an ergonomic and easy-to-use aid for sit-to-stand and transfers with turning. TurnSafe has been developed in close cooperation with medical personnel and users and therefore it efficiently looks out for the interests of both parties.

TurnSafe features a large foot plate with an anti-slip surface, adjustable height and lower leg supprts which are adjustable both sideways and in height. These features, together with the stable and grip-friendly handle, enable the user to actively participate in the sit-to-stand and transfer.

TurnSafe is relatively small and convenient, easy to have close at hand and very easy to rotate and manoeuvre, which makes it useful in most spaces and situations. By keeping one foot on the brake, the carer can avoid any unwanted movements while the user is rising up, standing or sitting down. The wheels, and the possibility to easily disassemble TurnSafe into three pieces, facilitates both transportation and storage.



## **TECHNICAL INFORMATION**

Material:	Powdercoated steel
Foot plate diameter:	Ø 390 mm/15,6"
Height of lower handle, for the user:	790-1130 mm/31,6"-42,2"
Height of upper handle, for the caregiver:	860-1200 mm/34,4"-48"
Height of lower leg supports:220-530 mm/8,8"-21,2", adjustable in steps of 25 mm/1"	
Weight:	16,5 kg/36,3 lbs
Safe working load:	150 kg /330 lbs

Product Art.No.
TurnSafe 400025544





## **ACCESSORIES**

Art.No.	Product
7215	Padding for leg support, 1 pair.

SystemRoMedic $^{7M}$  is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

 $\textbf{Transfer}, \ \text{between two locations}. \ \textbf{Positioning}, \ \text{within the same location}.$ 

Support, for mobility. Lifting, both manual and mechanical.



Handicare AB
Box 640, SE 175 27 Järfälla, SWEDEN
Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
www.handicare.com
info@handicare.se