


Patient's name: \_\_\_\_\_

Social insurance number: \_\_\_\_\_



StandingVest is specially designed for lifting with a standing shell. For the caregiver, it is easy to apply and it gives the patient a feeling of security. It is designed to provide upper-body support and is held in place by leg harness at the groin. StandingVest is designed to activate the patient safely and comfortably.



### Functional inspection

#### Visual and mechanical inspection

Check the condition and function of the sling regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. Apply load to the device and check to ensure that clasps, handles, etc. withstand heavy load. If there are signs of wear, the product must be discarded.

#### Which size is used?

Small       Medium       Large       X-Large

#### Which material is used?

Polyester       Other

#### Which loop alternatives are used?

\_\_\_\_\_

## Instructions for use

Can the vest be left behind the patient in the wheelchair?

Yes  No

Can the vest be used for transfers to/from the toilet?

Ja  Nej

## Have staff that use the assistive device:

Read the instructions?

Yes  No

Received practical instruction in patient transfers?

Yes  No

## Instructions to personnel:

e.g.; Placement of head/ arms/hands? Placement of legs/feet? Use bed functions.

---

---

---

---

Number of caregivers: \_\_\_\_\_

## Instructions to patient:

e.g.; In what way can the patient participate actively in the transfer situation? Which instructions are used?

---

---

---

---

## Other:

---

---

## Trial-fitted by:

Date: \_\_\_\_\_

Phone: \_\_\_\_\_

Phone: \_\_\_\_\_

Handicare AB  
Veddestav. 15, Box 640  
SE-175 27 Järfälla  
SWEDEN

Tel: +46 (0)8-557 62 200  
Fax: +46 (0)8-557 62 299  
E-mail: [info@handicare.se](mailto:info@handicare.se)

Internet: [www.handicare.com](http://www.handicare.com)

## In the events of changes or questions, contact the trial-fitter.

Manuals and instructions for functional testing are  
available on Handicare's website [www.handicare.com](http://www.handicare.com)

