

# MiniLift125

SystemRoMedic™

**SWL**  
125 kg  
275 lbs

*Easily adjustable lift arm*

*Easy to manoeuvre*

*Grip-friendly handles and side supports provide multiple hand placement options for the user*

*Quick guide on the lift*

**Unique sit-to-stand pattern**

*High capacity battery with long service life*

*Compact format*

*Built-in charging*

*Soft, adjustable lower leg support*

*Lockable rear castors*

*Footplate with low step-in height and anti-slip surface*



*Electrical base widening*

## A unique solution for active sit-to-stand – with electrical base widening

MiniLift125 is a mobile sit-to-stand lift which is used, in combination with a specially adapted lifting vest, to assist users with reduced function when rising from a sitting to a standing position. MiniLift125 features a construction offering a totally unique and extremely comfortable sit-to-stand experience. The lift is also compact and convenient and very easy to manoeuvre; features which, in combination with the possibility to adjust the settings for an optimal fit for every user, makes it very useful in most locations and situations.



# MiniLift125

## Extremely comfortable, active sit-to-stand with natural movement pattern

MiniLift125 has been developed to, as gently and comfortably as possible, assist the user when rising from a sitting to a standing position. The handles and the side supports offer multiple hand placement options for the user. The stable anti-slip foot plate is very low and the lower leg supports, as well as the lifting arm, can be adjusted to fit every user optimally. When MiniLift125 is used in combination with the appropriate lifting accessories, the user is supported under the feet, for the front of the lower legs and behind the back, which provides for a safe and secure, yet active sit-to-stand procedure. The construction safely moves the user forwards and upwards in a natural movement pattern and, at the same time, the leg muscles and balance are exercised. MiniLift125 is manufactured in steel for stability and strength and is suitable for users weighing up to 125 kg/275 lbs.

Active sit-to-stand with natural movement pattern.



## PRODUCT SPECIFICATION

MiniLift125 \_\_\_ Mobile sit-to-stand lift with electrical base widening \_\_\_ Art. No. 400641434

## TECHNICAL INFORMATION

Safe working load (SWL) \_\_\_\_\_ 125 kg / 275 lbs  
 User length, min-max \_\_\_\_\_ 1400-2000 mm / 55.1-78.7"  
 Material \_\_\_\_\_ Powder-coated steel  
 Height \_\_\_\_\_ 995 mm / 39.2"  
 Length \_\_\_\_\_ 900 mm / 35.4"  
 Width of base, outer measurement \_\_\_\_\_ 600-890 mm / 23.6-35.0"  
 Width of base, inner measurement \_\_\_\_\_ 445-745 mm / 17.5-29.3"  
 Weight \_\_\_\_\_ 39 kg / 86.0 lbs  
 Base widening \_\_\_\_\_ Electrical  
 Height of footplate, for user \_\_\_\_\_ 50 mm / 2.0"  
 Castor size, diameter \_\_\_\_\_ 80 mm / 3.1"  
 Lifting speed, without load \_\_\_\_\_ 35 mm / 1.4" per sec  
 Emergency lowering \_\_\_\_\_ Manual and electrical



MiniLift125 meets the requirements of the Medical Devices Directive for Class 1 products 93/42/EEC.



Ergonomic and easy-to-use hand control with 4 buttons.



New control box with diagnostics, e.g. lift counter and overload indicator.

Emergency stop button readily available on the control box.

## ACCESSORIES



Art. No.	Product name	Description
45500003	ThoraxSling	Polyester, X-Small
45500004	ThoraxSling	Polyester, Small
45500006	ThoraxSling	Polyester, Medium
45500007	ThoraxSling	Polyester, Large
45500008	ThoraxSling	Polyester, X-Large
45500009	ThoraxSling	Polyester, XX-Large
45550004	ThoraxSling	Wipeable, Small
45550006	ThoraxSling	Wipeable, Medium
45550007	ThoraxSling	Wipeable, Large
45550008	ThoraxSling	Wipeable, X-Large
45550009	ThoraxSling	Wipeable, XX-Large
45590004	ThoraxSling	Disposable, non-woven, Small, 5 pcs
45590006	ThoraxSling	Disposable, non-woven, Medium, 5 pcs
45590007	ThoraxSling	Disposable, non-woven, Large, 5 pcs
45590008	ThoraxSling	Disposable, non-woven, X-Large, 5 pcs
45600004	ThoraxSling	With seat support, polyester, Small
45600006	ThoraxSling	With seat support, polyester, Medium
45600007	ThoraxSling	With seat support, polyester, Large
45600008	ThoraxSling	With seat support, polyester, X-Large
45690004	ThoraxSling	With seat support, disposable, non-woven, Small, 5 pcs
45690006	ThoraxSling	With seat support, disposable, non-woven, Medium, 5 pcs
45690007	ThoraxSling	With seat support, disposable, non-woven, Large, 5 pcs
45690008	ThoraxSling	With seat support, disposable, non-woven, X-Large, 5 pcs
70200033	CalfStrap	For MiniLift125 as from serial number 125001653
770200089	Hand control HB33-6	6 buttons, with service, battery status and load indicators
70200029	External charger Linak	External charger for all models and variants of MiniLift

### ThoraxSling SystemRoMedic™

ThoraxSling is a specially adapted and very comfortable sit-to-stand vest that is used in combination with Handicare's SystemRoMedic™ MiniLift sit-to-stand lifts. It is easy for the caregiver to apply, while giving the user an added sense of security during rising and standing training. ThoraxSling also facilitates toileting. ThoraxSling is available in two models; one which provides support behind the back (ThoraxSling) and one that also gives support under the seat (ThoraxSling, with seat support), and in several materials and sizes.

### Other models and variants of MiniLift mobile sit-to-stand lifts

MiniLift is available in several models and variants with different lifting capacity, with electrical or manual base-width adjustment and with standard height or low legs for use with lower beds.

Art. No.	Product name	Description
401100334	MiniLift200	Mobile sit-to-stand lift with electrical base widening, SWL 200 kg/440 lbs
60300010	MiniLift160EM	Mobile sit-to-stand lift with manual base widening, SWL 160 kg/350 lbs
60300011	MiniLift160EM	Mobile sit-to-stand lift with manual base widening, low legs, SWL 160 kg/350 lbs
60300012	MiniLift160EE	Mobile sit-to-stand lift with electrical base widening, SWL 160 kg/350 lbs
60300013	MiniLift160EE	Mobile sit-to-stand lift with electrical base widening, low legs, SWL 160 kg/350 lbs

Please visit [www.handicare.com/us](http://www.handicare.com/us) for more information.



**handicare**

2201 Hangar Place, Suite 200  
 Allentown, PA 18109  
 Tel: 610 266 5260 • Fax: 610 266 5266  
 Toll Free: 866 276 LIFT/5438  
[www.handicare.com/us](http://www.handicare.com/us)  
[customerservice.usa@handicare.com](mailto:customerservice.usa@handicare.com)

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

**Transfer**, between two locations. **Positioning**, within the same location.

**Support**, for mobility. **Lifting**, both manual and mechanical.