

Technical drawings

Commode/shower chair soft seat


Article no:

LI2135.011 - Commode/shower chair soft seat, with closed front, 610 mm

LI2135.111 - Commode/shower chair soft seat, with closed front, 530 mm

LI2135.021 - Commode/shower chair soft seat, with open front, 610 mm

LI2135.121 - Commode/shower chair soft seat, with open front, 530 mm



	A	B	C
LI2135.011/LI2135.012	570 mm/22.4"	495 mm/19.5"	550 mm/21.7"
LI2135.111/LI2135.121	490 mm/19.3"	415 mm/16.3"	470 mm/18.5"

	mm	inch
D	555	21.9"
E	635	25"
F	450	17.7"
G	415	16.3"
H	490	19.3"
I	250	9.8"
J	625	24.6"